



Industrial Freezer Sales

A division of IDS

5311 Derry Avenue Building D Agoura Hills CA 91301

Phone (818) 597-4300 Fax (818) 597-4301

www.freezerlink.com

Security Guard Door Locks, Software



Security Guard is a keyless system. **There are no keys to lose or to share.** Security Guard comes with 250 user codes, allowing up to 250 users an individual access method - choose from Magstripe ID card, or Proximity ID card. Use existing identification cards. Security Guard **reads Magstripe Track 2 and HID® cards**, eliminating the need to buy new cards. Security Guard has a non-volatile memory, so stored information (users, supervisors, audit trails, etc.) is safe. Adding and deleting users is done through the unit **or optional Security Guard software.**

Secure Access to Valuable Samples

Secure Access to Valuable Samples Stored in Refrigerators and Freezers for Laboratory Scientific and Medical Storage Applications. Universal Design **may be applied to Glass Door or Solid Door** Units. Available in Your Choice of Personal Identification (PIN) Number Dual Systems: Proximity Reader/PIN/Keypad (**Reads HID Cards**) or Magstripe/ PIN/Keypad (**Reads Magstripe 2 Cards**).



Secure Access Features

- **250 User/Supervisor Codes**
- Audit Trail of Last 1500 Access Attempts
- Date, Time, and User Name.
- **Keyless Motor Drive Latching Mechanism,**
- Battery Operated with 4-AA Batteries (not Included)
- Low Battery Indicator.
- Jump Start 9-Volt Battery Power Port Allows the User to Access When the Batteries are expired.
- Audible Alarm Bolt System.
- **Non Volatile Memory.**
- **Green LED Light Confirms User Access.**
- Standard RJ11 Communications Port for Audit Trail Download.
- **Meets JCAHO Standards for use in medical facilities.**



Industrial Freezer Sales

A division of IDS


5311 Derry Avenue Building D Agoura Hills CA 91301

Phone (818) 597-4300 Fax (818) 597-4301

www.freezerlink.com

Access Rights Window

Allows the operator to choose which user can have access to each lock in the system. All locks and users in the system will appear. Can be re-arranged view by users showing which locks the user has access to and which locks the user does not have access to. Closed loop system of check unchecked boxes - indicates which entries on the access rights screen have been uploaded into the respective locks.



to

Lock Editor

Allows the operator to enter EAC-1000B serial numbers and pass codes. Enter the type of access hardware provided on the lock. Enter the amount of time that the operator wishes the lock to remain open.

IDS Part

Security Locks Pre-Installed and Field Kits

PKKPGPFI

Installed for Proximity Key/Key Pad (GP & UC Models)

MKKPGPFI

Installed for Magnetic Key/Key Pad (GP & UC Models)

PKKPPRFI

Installed for Prox Key/Key Pad (Prem, Non-GP Models) Door

MKKPPRFI

Installed for Magnetic Key/Key Pad (Prem., Non-GP Models) Door

PKKPGPKIT

Field for Prox Key/Key Pad (GP & UC Models)

MKKPGPKIT

Field for Mag Key/Key Pad (GP & UC Models)

PKKPPRKIT

Field for Prox Key/Key Pad (Prem., Non-GP Models) Per Door

MKKPPRKIT

Field for Mag Key/Key Pad (Prem., Non-GP Models) Per Door

LockSoft

Security Lock Software

As always, the best source of information for Industrial, Scientific and Medical Freezers and Refrigerators, equipment and practical, cost effective solutions is our 100 page website; www.freezerlink.com